

# The Alere Cholestech LDX® System Ordering Guide



## Test Cassettes (Box of 10)

Lipid Profile•GLU	10-991
Lipid Profile	10-989
TC•HDL•GLU	10-990
TC•HDL	10-987
TC•GLU	10-988
TC	10-986

Important note: Refrigerate upon arrival



## Controls/Calibration Verification

Multianalyte Controls Level 1 & 2 (1 x 2 Levels, 0.25mL each)	88773
Multianalyte Controls Level 1 & 2 (2 x 2 Levels, 2mL each)	88769
Cal Ver Material Lipids/GLU Levels 1-4 (1 x 4 Levels)	88770

Important note: Refrigerate upon arrival



## MiniPet™ Pipettes/Pipette Tips

40 µL Lipids/Glucose: All Samples	13-014
Pipette tips (50)	11-010



## Carrying Case

Carrying Case	10-973
---------------	--------



## Alere Universal Printer

Alere Universal Printer	14-973
Printer Cable Kit	11-787



## Capillary Tubes/Plungers

Capillary Tubes, 40 µL (50) (Lipids/Glucose)	52193
Capillary Plungers (50) (fits all capillary tubes)	10-311



## Lancets

Medlance® (200)	13-791
-----------------	--------



## White Printer Labels

For Alere Universal Printer (1 roll of 400)	11-785
---	--------



## Optics Check

Optics Check Cassette	10-228
-----------------------	--------

**Alere**™ Call 1.877.441.7440 or visit [www.alere.com](http://www.alere.com) today.